

## DinoLite support on MAC

**System requirement:**

PowerPC G3/G4/G5 or Intel-based Macintosh

**DinoLite model:**

DinoLite : AM211 / AM311x

DinoLite Plus: AM313x

DinoLite Pro: AM413x

DinoEye: AM423

		Minor Version									
MAC OS		.0	.1	.2	.3	.4	.5	.6	.7	.8	.9
Major Version	X10.0	No	No	No	No	No					
	X10.1	Ok	Ok	Ok	Ok	Ok	Ok				
	X10.2	Ok	Ok	Ok	Ok	Ok	Ok	Ok	Ok	Ok	
	X10.3	Ok*	Ok*	Ok*	Ok*	Ok*	Ok*	Ok*	Ok*	Ok*	
	X10.4	Ok	Ok	Ok	Ok	Ok	Ok	Ok	Ok	Ok	Ok
	X9	No	Ok	Ok							

(\*) QuickTime 6.5.1 or later is required.