

Dino-Lite AD4013MZTL

The AD series have a detachable cap and is delivered with both a closed and an open exchangeable cap. This MZTL model has a long working distance, polarisation filter and a third (exchangeable) polarisation cap as well as aluminum alloy housing. The enhanced color sensor provides superior color rendering and improved low light imaging.



Technical data

Model	AD4013MZTL (Long Working Distance)
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	10x ~ 92X
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√
LED color	White
Material	Full aluminum alloy body
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√ (line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	Polarisation
Compatibility	Windows XP, Vista, 7 and recent MAC OS
Software (included)	Windows: DinoCapture, Mac: DinoXcope
Other	best model for repair, rework or assembly tasks on reflective surfaces