

Dino-Lite AD413MT

The AD series have a detachable cap and is delivered with both a closed and an open exchangeable cap. The MT model additionally has an aluminum alloy housing.



Technical data

Model	AD413MT
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	10x~70x, 200x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√
LED color	White
Material	Full aluminum alloy body
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√(line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and recent MAC OS
Software (Included)	Windows: DinoCapture, Mac: DinoXcope
Other	designed for the professional market. Equipped with 2 exchangeable caps (open cap and closed cap)