

## Product data sheet

[www.dino-lite.eu](http://www.dino-lite.eu)

### Dino-Lite AD413T

The Dino-Lite AD series is designed for the professional market. The AD413T comes with 2 exchangeable caps, an open cap for most common applications and a closed cap designed to protect the lens and built-in LED lights from dust and moisture. Other caps for the AD series are available.



### Technical data

Model	AD413T
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	10x~70x, 200x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√
LED color	White
Material	Composite
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√(line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and recent MAC OS
Software (Included)	Windows: DinoCapture, Mac: DinoXcope
Other	designed for the professional market. Equipped with 2 exchangeable caps (open cap and closed cap)