

## Dino-Lite AD413TL

The Dino-Lite AD series is designed for the professional market. The AD413TL comes with 2 exchangeable caps, an open cap for most common applications and a closed cap designed to protect the lens and built-in LED lights from dust and moisture. Also available are the optional FC-L-MA1 mirror angle and FC-L-TB1 tube cap. This TL model has a long working distance making it an ideal solution for tasks where you need more working distance to the object, for example repair, rework or assembly tasks.



### Technical data

Model	AD413TL (Long Working Distance)
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	10x ~ 92X
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√
LED color	White
Material	Composite
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√ (line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and recent MAC OS
Software (included)	Windows: DinoCapture, Mac: DinoXcope
Special	best model for repair, rework or assembly tasks
Other	Optional: FC-L-MA1 mirror angle and FC-L-TB1 tube