

Dino-Lite AD413ZT

By using a special polarization filter, the AD413ZT is the right choice when working with high glare materials such as plastics and metals. This model is delivered with three exchangeable caps: an open cap, a closed cap and a cap with polarization filter.



Technical data

Model	AD413ZT
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	10x~70x, 200x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√
LED color	White
Material	Composite
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√ (line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	built-in polarizer
Compatibility	Windows XP, Vista, 7 and recent MAC OS
Software (Included)	Windows: DinoCapture, Mac: DinoXcope
Other	The polarization effect can be turned on/off or can be adjusted to offer full or half polarization effect.