

Dino-Lite AM211

The Dino-Lite AM-211 is the basic model USB microscope with two ranges of magnification, 10~70 and up to 200 x magnification.



Technical data

Model	AM211
Resolution	640x480pixels (VGA)
Magnification	10x~70x, 200x
Interface	USB 1.0
Nr. of LED's	4
LED on/off switchable	-
LED color	White
Material	Composite
Microtouch	-
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second
Measurement	-
Calibration	-
Filter/Diffuser	-
Compatibility	Windows: XP, Vista and 7, MAC: Driver only
Software (Included)	Windows: DinoCapture, MAC: Driver only
Special	Affordable solution for home or school use. Suitable for hobbies & collections (coins, stamps, watches, jewels, etc.)