

## Dino-Lite AM311S

The AM311S is the USB 2.0 version of the AM211. This model has 8 LED's and has the possibility to turn the LED's on and off in the DinoCapture application.



### Technical data

<b>Model</b>	<b>AM311S</b>
Resolution	640x480pixels (VGA)
Magnification	10x~70x, 200x
Interface	High speed USB image transmission, USB 2.0
Nr. of LED's	8
LED on/off switchable	√
LED color	White
Material	Composite
Microtouch	-
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	-
Calibration	-
Filter/Diffuser	-
Compatibility	Windows XP, Vista and 7. MAC: Driver only
Software (Included)	Windows: DinoCapture, MAC: Driver only
Special	Affordable solution for home or school use. Suitable for hobbies & collections (coins, stamps, watches, jewels, etc.)