

Dino-Lite AM313T

The AM313T has both the measurement functions and the MicroTouch button to capture pictures directly from the Dino-Lite.



Technical data

Model	AM313T
Resolution	640x480pixels (VGA)
Magnification	10x~70x, 200x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√
LED color	White
Material	Composite
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√(line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and recent MAC OS
Software (Included)	Windows: DinoCapture, Mac: DinoXcope
Special	Affordable solution for industrial, scientific and business use