

Dino-Lite AM313T5

The AM313T5 has a fixed 500x magnification for high power applications. To limit glare onto the object, the AM313T5 has diffused lighting.



Technical data

Model	AM313T5
Resolution	640x480pixels (VGA)
Magnification	500x fixed focus magnification
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√
LED color	White
Material	Composite
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√(line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and recent MAC OS
Software (Included)	Windows: DinoCapture, Mac: DinoXScope
Special	For examining very small objects when high magnification is needed. Scientific, biological or specialized industrial use.