

## Dino-Lite AM4013MT

Dino-Lite X Pro with aluminum alloy housing (with Microtouch). The aluminum housing offers enhanced protection and durability and has the best look and feel.



### Technical data

<b>Model</b>	<b>AM4013MT</b>
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	10x~70x, 200x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√
LED color	White
Material	Full aluminum alloy body
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√(line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and recent MAC OS
Software (Included)	Windows: DinoCapture, Mac: DinoXcope