

## Dino-Lite AM4013MZT

By using a special polarization filter, the AM4013MZT is the right choice when working with high glare materials such as plastics and metals. This model also sports an aluminum metal housing for enhanced protection and durability and the best look and feel.



### Technical data

<b>Model</b>	<b>AM4013MZT</b>
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	10x~70x, 200x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√
LED color	White
Material	Full aluminum alloy body
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√ (line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	built-in polarizer
Compatibility	Windows XP, Vista, 7 and recent MAC OS
Software (Included)	Windows: DinoCapture, Mac: DinoXcope
Special	the polarization effect can be switched on/off or can be adjusted to offer full or half polarization effect.