

Dino-Lite AM4013TL-M40

The AM413TL-M40 is a special version of the long working distance series with a lower magnification rate of between 5 and 40 times.



Technical data

Model	AM4013TL-M40
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	5~40x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√
LED color	White
Material	Composite
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√
Calibration	√
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and recent MAC OS
Software (included)	Windows: DinoCapture, Mac: DinoXcope