

Dino-Lite AM412MP

Aluminum alloy version of the Dino-Lite for television. The AM412MP has a specially designed nozzle with built in light diffuser to reduce glare onto the object.



Technical data

Model	AM412MP
Resolution	628x586 pixels
Magnification	10x~70x, 200x
Interface	TV PAL
Nr. of LED's	8
LED on/off switchable	-
LED color	White
Material	Full aluminium alloy body
Microtouch	-
Output	-
Video Frame Rate	25 frames-per-second
Measurement	-
Calibration	-
Filter/Diffuser	built in light diffuser
Compatibility	-
Software (Included)	-
Special	-
Other	-