

Dino-Lite AM412MZP

Aluminum alloy version for television with polarizer filter.



Technical data

Model	AM412MZP
Resolution	628x586 pixels
Magnification	10x~70x, 200x
Interface	TV PAL
Nr. of LED's	8
LED on/off switchable	-
LED color	White
Material	Full aluminium alloy body
Microtouch	-
Output	-
Video Frame Rate	25 frames-per-second
Measurement	-
Calibration	-
Filter/Diffuser	built-in polarizer
Compatibility	-
Software (Included)	-
Special	
Other	the polarization effect can be switched on/off or can be adjusted to offer full or half polarization effect.