

## Dino-Lite AM412PTL

Television version with long working distance and LED on-off.



### Technical data

<b>Model</b>	<b>AM412PTL</b>
Resolution	628x586 pixels
Magnification	10x~92X
Interface	TV PAL
Nr. of LED's	8
LED on/off switchable	√ (with microtouch button)
LED color	White
Material	Composite
Microtouch	-
Output	-
Video Frame Rate	25 frames-per-second
Measurement	-
Calibration	-
Filter/Diffuser	built-in light Source Diffuser
Compatibility	-
Software (Included)	-
Special	
Other	<p>Working Distance table:            92X Magnification -&gt; 4.7 cm/1.9" from object            79X Magnification -&gt; 5.0 cm/2.0" from object            52X Magnification -&gt; 6.3 cm/2.5" from object            27X Magnification -&gt; 11.5 cm/4.5" from object            20X Magnification -&gt; 14.4 cm/5.7" from object</p>