

## Dino-Lite AM413FVT2

The AM413FVT2 has special lighting with LED's emitting UV light around the 365 nm spectrum.



### Technical data

<b>Model</b>	<b>AM413FVT2</b>
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	10x~70x, 200x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	4
LED on/off switchable	√
LED color	4 UV LED's 365nm (Emission Filter Cutoff Wavelength: 425nm)
Material	Composite
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√(line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and recent MAC OS
Software (Included)	Windows: DinoCapture, Mac: DinoXcope
Special	For scientific and engineering research/inspection. AM413FVT2 is suggested for Law enforcement