

Dino-Lite AM413MT-VW

The VW series has two types of LED lights that are *switchable* between white and UV light. The Dino-X-Lite models have the added advantage of an aluminum metal housing which offers enhanced protection and durability and has the best look and feel.



Technical data

Model	AM413MT-VW
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	10x~70x, 200x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	4+4
LED on/off switchable	√
LED color	4 UV LED's 390nm-400 nm + 4 White LED's (Emission Filter Cutoff Wavelength: 425nm)
Material	Full aluminum alloy body
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√ (line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and recent MAC OS
Software (Included)	Windows: DinoCapture, Mac: DinoXcope
Special	For scientific and engineering research/inspection