

Dino-Lite AM413MT5

The AM413MT5 has a aluminum alloy body and features fixed 500x magnification for high power applications. To limit glare onto the object, the AM413MT5 has diffused lighting. The aluminum housing offers enhanced protection and durability and has the best look and feel.



Technical data

Model	AM413MT5
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	500x fixed focus magnification
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√
LED color	White
Material	Full aluminum alloy body
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√ (line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	Light Source Diffuser
Compatibility	Windows XP, Vista, 7 and recent MAC OS
Software (included)	Windows: DinoCapture, Mac: DinoXcope
Special	For examining very small objects when high magnification is needed. Scientific, biological or specialized industrial use.