

Dino-Lite AM413MTL-FVW

The VW series has two types of LED lights that are switchable between white and UV light. This model has a long working distance and a filter to filter out the returning UV light so that an image with natural colors can be obtained. The aluminum housing offers enhanced protection and durability and has the best look and feel.



Technical data

Model	AM413MTL-FVW
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	10x ~ 92X
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	4+4
LED on/off switchable	√
LED color	4 UV LED's 390nm-400 nm + 4 White LED's (Emission Filter Cutoff Wavelength: 425nm)
Material	Full aluminum alloy body
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√ (line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	√, UV filter
Compatibility	Windows XP, Vista, 7 and recent MAC OS
Software (Included)	Windows: DinoCapture, Mac: DinoXcope
Special	For scientific and engineering research/inspection.