

## Dino-Lite AM413MTL-VW

The VW series has two types of LED lights that are *switchable* between white and UV light. The Dino-X-Lite models have the added advantage of an aluminum metal housing which offers enhanced protection and durability and has the best look and feel. This model also has a long working distance.



### Technical data

| Model                 | AM413MTL-VW   |
|-----------------------|---|
| Resolution            | 1,3 megapixel (1280x1024 pixels)  |
| Magnification         | 10x ~ 92X   |
| Interface             | USB 2.0, High speed USB transmission  |
| Nr. of LED's          | 4+4   |
| LED on/off switchable | √   |
| LED color             | 4 UV LED's 390nm-400 nm + 4 White LED's<br>(Emission Filter Cutoff Wavelength: 425nm) |
| Material              | Full aluminum alloy body  |
| Microtouch            | √   |
| Output                | Picture, video and time-lapsed video  |
| Video Frame Rate      | up to 30 frames-per-second, 15fps at 1.3M   |
| Measurement           | √ (line, radius, circle, 3-point circle, etc.)  |
| Calibration           | √   |
| Filter/Diffuser       | -   |
| Compatibility         | Windows XP, Vista, 7 and recent MAC OS  |
| Software (Included)   | Windows: DinoCapture, Mac: DinoXcope  |
| Special               | For scientific and engineering research/inspection.                                   |