

Dino-Lite AM413MTL

Dino-Lite X Pro with aluminum alloy housing and long working distance and Microtouch. The aluminum housing offers enhanced protection and durability and has the best look and feel.



Technical data

Model	AM413MTL
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	10x ~ 92X
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√
LED color	White
Material	Full aluminum alloy body
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√ (line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and recent MAC OS
Software (included)	Windows: DinoCapture, Mac: DinoXcope
Special	best model for repair, rework or assembly tasks