

## Dino-Lite AM413T

For the serious precision applications the Dino-Lite Pro AM413T has a high resolution of 1,3 megapixel and an adjustable magnification of appr. 10 to 70 times as well as 200x. This model also provides the possibility to calibrate the microscope for measurements. Most used model in the industrial field.



### Technical data

<b>Model</b>	<b>AM413T</b>
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	10x~70x, 200x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√
LED color	White
Material	Composite
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√(line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and recent MAC OS
Software (Included)	Windows: DinoCapture, Mac: DinoXcope