

Dino-Eye AM422PX

The Dino-Eye is a range of digital eyepiece cameras. You can replace the ocular (eyepiece) of regular microscopes with the Dino-Eye which turns the regular microscope into a digital microscope. The AM422PX with TV output can be directly connected to any TV monitor and has different adapters to fit microscopes with 23, 30 or 30.5 mm eyepiece.



Technical data

Model	AM422PX Dino-Eye TV
Resolution	510×496 pixels (PAL)
Magnification	72,5 times the magnification of objective
Interface	TV Output (PAL)
Nr. of LED's	-
LED on/off switchable	-
LED color	-
Material	Aluminium
Microtouch	-
Output	-
Video Frame Rate	up to 25 frames-per-seconds
Measurement	-
Calibration	-
Filter/Diffuser	-
Compatibility	-
Software (Included)	-
Special	The AM422PX has different adapters to fit microscopes with 23, 30 or 30,5 mm eyepieces