


## Dino-Lite MS15X X/Y table

The MS15X is an advanced XY table that can be used independently or in combination with the MS35B or MS36B column stands.

The MS15X X/Y table base with detachable rotating table is designed for use with all sorts of applications that require fine examination in an X- and Y-axis movement. With the 360 degree rotating table, one can easily carry out all kinds of precision examinations.

The MS15X can be used stand-alone or attached to the MS35B and MS36B stands (with the included C-Clamps)



### Technical data

<b>Model</b>	<b>MS15X</b>
Unit dimensions (WxLxH)	16,5 x 13,5 x 7,5 cm
Unit weight	appr. 1 kg
Other	Can be used independently or mounted on MS35B/MS36B stand