

## Dino-Lite MS52B stand

Dino-Lite heavy duty jointed flex arm stand.

The MS52B is a jointed articulating flex-arm type stand with holster for use with the Dino-Lite hand held digital microscopes. The precision machine crafted components provide smooth adjustment and unlimited horizontal positioning in a well balanced ergonomic design. The MS52B can be easily clamped onto and secured to any workstation or table top. Manufactured from durable, lightweight recycled aluminum and able to be folded flat makes the versatile MS52B stand conveniently portable.

Note: The available MS52BA2 option includes the MS52B flex-arm stand with an additional rigid vertical floor stand and base assembly for use when a workstation table top is not convenient or available.



### Technical data

Model	MS52B
Unit dimensions (WxLxH)	36 x 15,3 x 5 cm
Unit weight	appr. 1,5 kg
Horizontal & vertical work area	621 - 921 mm
Top swivel	360 degrees
Other	Optional vertical base MS52BA2 is available