

## Dino-Lite MS52BA2 stand

Heavy duty table base for the MS52B.

The MS52BA2 is a combination of the rigid MS52B articulating flex-arm and a desktop base rack that holds the MS52B on top of a table. The desktop base rack is designed for situations where the standard clamp cannot be used to secure the MS52B on a stable surface.

The flex-arm is moved and adjusted with tightening knobs and clamps for positioning the arm to the area you wish to examine. Two of the flex-arm axis can be rotated 360 degrees and extended for further reach.

The MS52BA2 flex-arm is manufactured from durable, lightweight recycled aluminum and has the ability to be folded flat, which makes the MS52BA2 stand portable and ready for travel.



### Technical data

Model	MS52BA2
Unit dimensions (WxLxH)	36 x 15,3 x 5 cm
Unit weight	appr. 2270 g
Horizontal & vertical work area	621 - 921 mm
Top swivel	360 degrees
Other	-