



Compatibility Chart Mac OS X with DinoXcope

Dino-Lite Industrial

 The Industry Standard	Dino-Lite 211/311S/313/411/413/423 series	Dino-Lite 2xxx/3xxx series	Dino-Lite 3613/3713/3715 series	Dino-Lite 4013/4113/4023 series	Dino-Lite 4115/4515/4815 series (EDR/EDOF is unavailable)	Dino-Lite 7013/7023 series	Dino-Lite 7115/7515/7815**/7915**/7025X (EDR/EDOF is unavailable)	Dino-Lite 73115/73515/73915** series (EDR/EDOF is unavailable)
Mac OS X 10.5 (Leopard)	Version 1.16 + Driver	Version 1.16	no	Version 1.16	no	Version 1.16	no	no
Mac OS X 10.6 (Snow Leopard)	Version 1.16 + Driver	Version 1.16	no	Version 1.16	no	Version 1.16	no	no
Mac OS X 10.7 (Lion)	Version 1.16 + Driver	Version 1.16	no	Version 1.16	no	Version 1.16	no	no
Mac OS X 10.8 (Mountain Lion)	Version 1.16 + Driver	Version 1.16	no	Version 1.16	no	Version 1.16	no	no
Mac OS X 10.9 (Mavericks)	Version 1.19.2 + Driver	Version 1.19.2	no	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2
Mac OS X 10.10 (Yosemite)	Version 1.19.2 + Driver	Version 1.19.2	no	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2
Mac OS X 10.11 (El Capitan)	no	Version 1.20.4	no	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2
Mac OS X 10.12 (Sierra)	no	Version 1.20.4	no	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4
Mac OS X 10.13 (High Sierra)	no	Version 1.20.4	no	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4
Mac OS X 10.14 (Mojave)	no	Version 1.20.4	no	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4

Dino-Lite Medical

	Dino-Lite Medical DermaScopes MEDL4DW, MEDL7DW	Dino-Lite Medical DermaScopes MEDL4DM, MEDL7DM	Dino-Lite Medical EarScopes MEDL3E, MEDL4E, MEDL4EP	Dino-Lite Medical TrichoScopes, MEDL3H, MEDL4HM, MEDL4HMA, MEDL4HMA5, MEDL4HVW, MEDL7HM, MEDL7HMA, MEDL7HM4	Dino-Lite Medical CapillaryScopes, MEDL4N Pro*, MEDL4N5 Pro*	Dino-Lite Medical CapillaryScope MEDL4N5	Dino-Lite Medical IriScope MEDL4R	Dino-Lite Medical MEDL4FW
Mac OS X 10.5 (Leopard)	no	Version 1.16	Version 1.16	Version 1.16	no	Version 1.16	Version 1.16	no
Mac OS X 10.6 (Snow Leopard)	no	Version 1.16	Version 1.16	Version 1.16	no	Version 1.16	Version 1.16	no
Mac OS X 10.7 (Lion)	no	Version 1.16	Version 1.16	Version 1.16	no	Version 1.16	Version 1.16	no
Mac OS X 10.8 (Mountain Lion)	no	Version 1.16	Version 1.16	Version 1.16	no	Version 1.16	Version 1.16	no
Mac OS X 10.9 (Mavericks)	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2
Mac OS X 10.10 (Yosemite)	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2
Mac OS X 10.11 (El Capitan)	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2	Version 1.19.2
Mac OS X 10.12 (Sierra)	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4
Mac OS X 10.13 (High Sierra)	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4
Mac OS X 10.14 (Mojave)	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4	Version 1.20.4

\* 4515/7115/7515 models: FLC and AMR compatible from version 1.20.4.

\*\* 4815/7815/7915/73915 models: FLC and AMR compatible from version 1.20.4. No EDOF/EDR function.