

# MS08B



**Brand:** Dino-Lite

**Product Code:** MS08B

## Short Description

Compact basic stand

## Description

The MS08B is an ultra-portable mini stand that holds the Dino-Lite in place whenever you need to. The MS08B is a great little accessory to enhance the Dino-Lite's microscopic portability solution.

## Specification

| Compatibility        |                                                                 |
|----------------------|-----------------------------------------------------------------|
| Compatible with:     | All Dino-Lite models, except Extra Long Working Distance models |
| Housing              |                                                                 |
| Dimensions           | 7.5cm (W) x 10cm (L) x 14,5cm (H)                               |
| Weight               | 500 g                                                           |
| Features             |                                                                 |
| ESD safe             | No                                                              |
| Information          |                                                                 |
| Package contents     | Stand, HD-P1 holder                                             |
| Warranty information | 2 years European warranty                                       |
| Accessory range:     | Basic stands                                                    |
| Price range          | €50,00 - €150,00                                                |

**Product Gallery**

