Dino-Lite Europe - Phone: - Email:

TrichoScope Basic (MEDL3H)

Código del producto: MEDL3H

Short Description

VGA (640x480) Resolution UV model for Trichoscopy Up to 200x magnification Medical Device Class 1

Descripción

The Dino-Lite TrichoScope Basic (MEDL3H) is a great instrument for the examination of the human scalp and for hair analysis, e.g. for hair loss research. The special cone shaped cap allows easy separation of hair on the human head. With the high magnification of up to 200 times a single hair can be viewed with great detail



USB 2.0







Especificación

Lighting		
Light/ LED type	White	
Number of LEDs	8	
LED on/off switchable:	Yes	
Infrared filter	IR cut-filter >650 nm	
Diffuser available	No	
Polarizer	No	
Optics		
Magnification	10-70x, 200x	
Lens type	Glass with anti-reflection coating	
Sensor		
Sensor type	CMOS	
Resolution	VGA (640x480)	

Maximum frame rate	30 fps
Compatibility	
Interface	USB 2.0
Operating system	Windows 7, 8, 10 & 11, MacOS 10.9 and up
Software	DinoCapture 2.0 (Windows), DinoXcope (Mac OS)
Supported image formats	BMP, GIF, PNG, JPG, TIF, RAS, PNM, TGA,
(Windows)	PCX, MNG, WBMP, JP2, JPC, PGX
Supported video formats	WMV, FLV, SWF
(Windows)	
Supported image formats	JPEG, PNG
(MacOS)	
Supported video formats	MOV
(MacOS)	
Imaging standards	DirectShow, UVC
Housing	
Housing material	Composite/plastic housing
Magnification lock	No
Dimensions	10.5cm (L) x 3.2cm (D)
Weight	95 g
Cable length	1.8m
Features	
Special feature	No
Measurement	No
Calibration	No
Microtouch sensor	Yes
Information	
Package contents	Microscope, carry pouch, software CD, ZT-Z-CC1
	cap closure, user manual
Warranty information	2 years European warranty
Regulatory approval	Medical Device Class 1 – Medical Devices
	Regulation (EU) 2017/745
Price range	€300,00 - €450,00

Product Gallery



