

# RK-06F



**Marque:** Dino-Lite  
**Code produit:** RK-06F

## **Courte description**

Portable rack  
Fast assembly  
Fine-focus adjustment

## **Description**

The Dino-Lite RK-06F stand is a sturdy and stable mid-range stand solution, with fine focus adjustment that can be used with all Dino-Lite digital microscopes. It is constructed of stainless steel and lightweight aluminum and offers precise fine-focus adjustment as well as a quick vertical release function.

Its adjustable horizontal arm has a 15cm range with 360-degree rotation and with the RK-10-EX arm (excluded) its effective range can be extended to observe larger objects. Various optional accessories are available for the RK-06F, such as the RK-10-FX flexible arm or the RK-10-PX, XY arm. With these features, the RK-06F is the ideal mid-range stand solution for precise orientation and maximum flexibility.

## [AVAILABLE EXTENSIONS FOR THE RK-06A STAND](#)


---

### **RK-10-EX**




additional horizontal arm (included with RK-10A)

**RK-10-FX**




flexible arm extension

**RK-10-VX**



vertical arm extension

**RK-10-PX**



*(Temporarily unavailable)*

### Spécification

Compatibility	
Compatible with:	All Dino-Lite models
Housing	
Dimensions	22 cm (L) x 15 cm (W) x 27 cm (H)
Weight	1060 g
Features	
ESD safe	No

Information	
Package contents	Stand, HD-P1 holder, RK-10-EX horizontal arm
Warranty information	2 years European warranty
Accessory range:	Professional stands
Price range	€250,00 - €400,00

## Galerie du produit

