

RK-04F



Marca: Dino-Lite
Codice del prodotto: RK-04F

Short Description

Portable rack
Fast assembly
Fine-focus adjustment

Descrizione

The Dino-Lite RK-04F stand has a robust, lightweight and compact rack design. Additionally it has an integrated fine focus adjustment knob to make focusing easier and more precise. The RK-04F is the perfect accessory for users requiring flexibility.

All Dino-Lite microscopes are compatible with this stand with the exception of the Long Working Distance microscopes.

The package includes the stand and a HD-P1 holder.

[Also available: RK-04 without fine focus adjustment](#)

Specifiche

| Compatibility | |
|------------------|---|
| Compatible with: | All Dino-Lite models, except Long Working Distance models |
| Housing | |

| | |
|----------------------|----------------------------------|
| Dimensions | 14cm (W) x 11cm (L) x 16.5cm (H) |
| Weight | 630 g |
| Features | |
| ESD safe | No |
| Information | |
| Package contents | Stand, HD-P1 holder |
| Warranty information | 2 years European warranty |
| Accessory range: | Professional stands |
| Price range | €150,00 - €250,00 |

Product Gallery

