


BL-CDW



Merk: Dino-Lite
Productcode: BL-CDW

Korte omschrijving

Brightfield/ Darkfield Backlight illuminator

Omschrijving

The Dino-Lite backlight stage BL-CDW allows two different types of lighting: brightfield and darkfield. In the brightfield mode, the light is transmitted from below the sample that shows up vividly showing the form and structure of the specimen on a bright background. With darkfield illumination, the transmitted light will not reach the microscope and will therefore not be part of the image.

The image will show bright objects on an almost black background. This type of light is not only suitable for prepared specimen but also for live specimen such as flagellum, the pseudopod movement, spirochete and more. Both types of illumination allow for high contrast, clear images and adds a whole new dimension to the Dino-Lite digital handheld microscopes with magnifications over 100x.

Specificatie

Compatibiliteit	
Compatibel met:	All Dino-Lite models, except long working distance models
Behuizing	
Afmetingen	5.1cm (H) x 5.43cm (D)

Gewicht	178 g
Kenmerken	
ESD veilig	No
Informatie	
Inhoud van de verpakking	BL-CDW, power adapter
Garantie-informatie	2 years European warranty
accessoires:	Light and control
Prijsklasse	€250,00 - €400,00

Product galerij

