

CA1070

Merk: Dino-Lite

Productcode: CA1070



Korte omschrijving

Sturdy watertight carrying case

Omschrijving

Sturdy, watertight case for Dino-Lite microscopes. This heavy duty carrying case with foam inlay protects your Dino-Lite against dust, scratches and moisture. Watertight according to IP67 standard.

Specificatie

Compatibiliteit	
Compatibel met:	All Dino-Lite models
Behuizing	
Afmetingen	Exterior: 210 x 167 x 90mm Interior: 186 x 123 x 75mm
Gewicht	570 g
Kenmerken	
ESD veilig	No
Informatie	
Inhoud van de verpakking	Case, foam inlay
Garantie-informatie	2 years European warranty
accessoires:	Various
Prijsklasse	€0 - €50,00

Product galerij

