

CS-41

Merk: Dino-Lite

Productcode: CS-41



Korte omschrijving

Glass calibration slide

Omschrijving

The CS-41 is a calibration slide that can be used to calibrate your Dino-Lite Edge digital microscope or to be used with the Dino-Eye Edge eyepiece cameras. The total length is 10mm (10.000um) and the smallest gradation on the slide is 0.1mm (100um). The CS-41 calibration slide also contains QR codes in 3 different sizes that can be used for the Auto Calibration procedure within DinoCapture 2.0 software version 1.5.25 or newer. Auto calibration provides an easy way to calibrate a Dino-Lite Edge serie microscope. The DinoCapture software identifies the special QR codes and performs the calibration automatically.

Row A: magnifications of 120x and above

Row B: magnifications of 40x - 220x

Row C: magnifications of 20x - 50x

[Click here to download the Auto Calibration procedure](#)

Specificatie

Compatibiliteit	
Compatibel met:	All Dino-Lite models

Behuizing	
Afmetingen	75mm (L) x 25mm (W) x 1,5mm (H)
Gewicht	10 g
Kenmerken	
ESD veilig	No
Informatie	
Inhoud van de verpakking	1 piece
Garantie-informatie	2 years European warranty
accessoires:	Various
Prijsklasse	€0 - €100,00