

AM4023U

Merk: Dino-Lite

Productcode: AM4023U



Korte omschrijving

1.3 Megapixel resolution

Fits over ocular up to 36 mm

Advanced software, including measurements

Omschrijving

The Dino-Eye is a range of digital eyepiece cameras. You can replace the ocular (eyepiece) of regular microscopes with the Dino-Eye which turns the regular microscope into a digital microscope. The AM4023U is an universal eyepiece version that fits over existing ocular up to 36 mm.



USB 2.0



Measurement
functionality



1,3 Megapixels
1280 x 1024



Fits over oculars
up to 36 mm

Specificatie

Optiek

Magnification	Yes, built-in
---------------	---------------

Type lens	Glass with Anti-Reflection coating
-----------	------------------------------------

Sensor

Sensortype	CMOS
Resolutie	1.3 Megapixel (1280x1024)
Maximale framesnelheid	30 fps
Compatibiliteit	
Interface	USB 2.0
Besturingssysteem	Windows 7, 8, 10 & 11, MacOS 10.9 and up
Software	DinoCapture 2.0 (Windows), DinoXcope (Mac OS)
Ondersteunde beeldformaten (Windows)	BMP, GIF, PNG, JPG, TIF, RAS, PNM, TGA, PCX, MNG, WBMP, JP2, JPC, PGX
Ondersteunde videoformaten (Windows)	WMV, FLV, SWF
Ondersteunde beeldformaten (MacOS)	JPEG, PNG
Ondersteunde videoformaten (MacOS)	MOV
Beeldvormingsnormen	DirectShow, UVC
Wifi	No
Behuizing	
Materiaal van de behuizing	Metal body
Afmetingen	70mm (L) x 45 mm (D)
Gewicht	190 g
Kabellengte	1.8 m
Ocular outer diameter:	Up to 36 mm
Mounting type:	Over ocular/tube
Kenmerken	
Meting	Yes
Kalibratie	Yes
Microtouch-sensor	No
Informatie	
Inhoud van de verpakking	DinoEye eyepiece camera, carry pouch, software CD, calibration target, user manual
Garantie-informatie	2 years European warranty
Regelgevende goedkeuring	CE, FCC, ROHS
accessoires:	DinoEye Eyepieces
Prijsklasse	€250,00 - €500,00

Product galerij

