

EarScope Basic Wireless (MEDLW3E)

Merk: Dino-Lite Medical
Productcode: MEDLW3E








Korte omschrijving

640x480 Resolution
 Basic model for otoscopy
 20-30x Magnification
 Wireless
 Medical Device Class 1

Omschrijving

The EarScope Basic Wireless (MEDLW3E) is a compact and safe solution to look at the ear. The cap is detachable for hygienical reasons. The MEDLW3E has a 640x480 resolution and magnification of 20 ~ 30x. This model is delivered with disposable specula of 3, 4 and 5 mm (10 each).

 <p>USB 2.0</p>	 <p>Adjustable ~20x-30x</p>	 <p>VGA resolution 640 x 480</p>
 <p>4 White LEDs non switchable</p>	 <p>Wireless</p>	

Specificatie

Verlichting	
Licht / LED-type	White

Aantal leds	4
LED aan / uit schakelbaar:	Yes
Infrarood filter	IR cut-filter >650 nm
Diffusor beschikbaar	No
Polarisator	No
Optiek	
Magnification	20-30x
Type lens	Glass with anti-reflection coating
Sensor	
Sensortype	CMOS
Resolutie	VGA (640x480)
Maximale framesnelheid	30 fps
Compatibiliteit	
Interface	USB 2.0
Besturingssysteem	Windows 7, 8, 10 & 11, MacOS 10.9 and up
Software	DinoCapture 2.0 (Windows), DinoXcope (Mac OS)
Ondersteunde beeldformaten (Windows)	BMP, GIF, PNG, JPG, TIF, RAS, PNM, TGA, PCX, MNG, WBMP, JP2, JPC, PGX
Ondersteunde videoformaten (Windows)	WMV, FLV, SWF
Ondersteunde beeldformaten (MacOS)	JPEG, PNG
Ondersteunde videoformaten (MacOS)	MOV
Beeldvormingsnormen	DirectShow, UVC
Wifi	Yes
Behuizing	
Materiaal van de behuizing	Composite/plastic housing
Vergrotingslot	No
Afmetingen	10.5cm (L) x 3.2cm (D)
Gewicht	95 g
Kabellengte	1.8m
Kenmerken	
Meting	No
Kalibratie	No
Microtouch-sensor	No
Informatie	
Inhoud van de verpakking	Microscope, carry pouch, 10x 3 mm specula, 10x 4 mm specula, 10x 5 mm specula, software CD, user manual
Garantie-informatie	2 years European warranty

Regelgevende goedkeuring	Medical Device Class 1 – Medical Devices Regulation (EU) 2017/745
Prijsklasse	€500,00 - €700,00

Product galerij

