

# RK-02

**Merk:** Dino-Lite

**Productcode:** RK-02



## Korte omschrijving

Gooseneck stand

## Omschrijving

The RK-02 is a durable gooseneck stand solution designed for flexible and versatile usage with a Dino-Lite digital microscope. The C-clamp of the RK-02 fits to virtually any surface of up to 7 cm (~3 inch) thick, and can be mounted horizontally or vertically. With this sturdy stand, you can easily position the Dino-Lite by bending the 60 cm long arm and adjusting the rotating head. The base of the RK-02 can also be rotated providing even more freedom and it can be fixed with the quick release fastener for use on uneven surfaces. The RK-02 gooseneck arm is the ideal choice to flexibly use any of the Dino-lite models, except the high magnification models.

## Specificatie

| Compatibiliteit |                                                        |
|-----------------|--------------------------------------------------------|
| Compatibel met: | All Dino-Lite models, except High Magnification models |
| Behuizing       |                                                        |
| Afmetingen      | Gooseneck 60cm (L) x 1.5cm (D)                         |
| Gewicht         | 1.5 kg                                                 |
| Kenmerken       |                                                        |
| ESD veilig      | No                                                     |

| Informatie               |                           |
|--------------------------|---------------------------|
| Inhoud van de verpakking | Stand, HD-P1 holder       |
| Garantie-informatie      | 2 years European warranty |
| accessoires:             | Professional stands       |
| Prijsklasse              | €100,00 - €200,00         |

## Product galerij

