

DermaScope® Polarizer (MEDL4DW)



Merk: Dino-Lite Medical

Productcode: MEDL4DW

Korte omschrijving

1.3 Megapixel resolution

Versatility due to exchangeable caps

10-50x Magnification

Polarizer

Medical Device Class 1

Omschrijving

The Dino-Lite DermaScope® Polarizer (MEDL4DW) has a 1.3 megapixel camera that can capture sharp images of the skin and features a built-in and fully adjustable polarization filter to greatly reduce the gloss effect of the skin.

The Dino-Lite DermaScope® Polarizer is suitable for any physician or dermatologist or skin specialist. For those who put even more demands on the equipment and want to have the latest features and capabilities, the DermaScope® Polarizer HR is also available (MEDL7DW)



USB 2.0



Exchangeable
cap



Adjustable
~10-50x



Magnification
lock



1.3 megapixel
1280 x 1024



8 White LEDs
switchable



Polarizer
anti-reflection



Macrozoom
FOV & DOF

Specificatie

Verlichting	
Licht / LED-type	White
Aantal leds	8
LED aan / uit schakelbaar:	Yes
Infrarood filter	IR cut-filter >650 nm
Diffusor beschikbaar	No
Polarisator	Yes, linear
Optiek	
Magnification	10-50x
Type lens	Glass with anti-reflection coating
Sensor	
Sensortype	CMOS
Resolutie	1.3 Megapixel (1280x1024)
Maximale framesnelheid	30 fps
Compatibiliteit	
Interface	USB 2.0
Besturingssysteem	Windows 7, 8, 10 & 11, MacOS 10.9 and up
Software	DinoCapture 2.0 (Windows), DinoXcope (Mac OS)
Ondersteunde beeldformaten (Windows)	BMP, GIF, PNG, JPG, TIF, RAS, PNM, TGA, PCX, MNG, WBMP, JP2, JPC, PGX
Ondersteunde videoformaten (Windows)	WMV, FLV, SWF
Ondersteunde beeldformaten (MacOS)	JPEG, PNG
Ondersteunde videoformaten (MacOS)	MOV
Beeldvormingsnormen	DirectShow, UVC
Behuizing	
Materiaal van de behuizing	Composite/plastic housing
Vergrotingsslot	Yes
Afmetingen	10.5 cm (L) x 3.2 cm (D)
Gewicht	105 g
Kabellengte	1.8m

Kenmerken	
Bijzonderheid	No
Meting	No
Kalibratie	No
Microtouch-sensor	Yes
Informatie	
Inhoud van de verpakking	Microscope, carry pouch, software CD, user manual, N3C-C- Closed cap, N3C-L- Long cap
Garantie-informatie	2 years European warranty
Regelgevende goedkeuring	Medical Device Class 1 – Medical Devices Regulation (EU) 2017/745
Prijsklasse	€500,00 - €700,00

Product galerij

