Dino-Lite Europe - Phone: - Email:

RK-05F



Product Code: RK-05F

Short Description

Portable rack
Tilt and stay feature (30 degrees)
Fast assembly
Fine-focus adjustment

Description

RK-05F is a light-weight, ultra-portable stand with a one-step assemble/disassemble design.

If you want to observe your object or specimen from a different perspective, the Dino-Lite RK-05F stand provides a multi-angle view with its 30 degrees inclinable pole feature. Additionally it has an integrated fine focus adjustment knob to make focusing easier and more precise.

RK-05F is a light-weight, ultra-portable stand with a one-step assemble/disassemble design.

All Dino-Lite microscopes are compatible with this stand with the exception of the Long Working Distance microscopes.

The package includes the stand and a HD-P1 holder.

Also available: RK-05 without fine focus adjustment

Specification

All Dino-Lite models, except Long Working
Distance models
14cm (W) x 11cm (L) x 16.5cm (H)
636 g
No
Stand, HD-P1 holder
2 years European warranty
Professional stands
€200,00 - €300,00

Product Gallery

