


# CS-41

**Brand:** Dino-Lite  
**Product Code:** CS-41



## Short Description

Glass calibration slide

## Description

The CS-41 is a calibration slide that can be used to calibrate your Dino-Lite Edge digital microscope or to be used with the Dino-Eye Edge eyepiece cameras. The total length is 10mm (10.000um) and the smallest gradation on the slide is 0.1mm (100um). The CS-41 calibration slide also contains QR codes in 3 different sizes that can be used for the Auto Calibration procedure within DinoCapture 2.0 software version 1.5.25 or newer. Auto calibration provides an easy way to calibrate a Dino-Lite Edge serie microscope. The DinoCapture software identifies the special QR codes and performs the calibration automatically.

Row A: magnifications of 120x and above

Row B: magnifications of 40x - 220x

Row C: magnifications of 20x - 50x

[Click here to download the Auto Calibration procedure](#)

## Specification

|                  |                      |
|------------------|----------------------|
| Compatibility    |                      |
| Compatible with: | All Dino-Lite models |

| Housing              |                                 |
|----------------------|---------------------------------|
| Dimensions           | 75mm (L) x 25mm (W) x 1,5mm (H) |
| Weight               | 10 g                            |
| Features             |                                 |
| ESD safe             | No                              |
| Information          |                                 |
| Package contents     | 1 piece                         |
| Warranty information | 2 years European warranty       |
| Accessory range:     | Various                         |
| Price range          | €0 - €100,00                    |