

CA1070

Marka: Dino-Lite
Ürün Kodu: CA1070



Short Description

Sturdy watertight carrying case

Açıklama

Sturdy, watertight case for Dino-Lite microscopes. This heavy duty carrying case with foam inlay protects your Dino-Lite against dust, scratches and moisture. Watertight according to IP67 standard.

Özellikler

Compatibility	
Compatible with:	All Dino-Lite models
Housing	
Dimensions	Exterior: 210 x 167 x 90mm Interior: 186 x 123 x 75mm
Weight	570 g
Features	
ESD safe	No
Information	
Package contents	Case, foam inlay
Warranty information	2 years European warranty
Accessory range:	Various
Price range	€0 - €50,00

Product Gallery

